

	CLASS	COMPOUND	METABOLINK™	CUTOFF (ng/mL)	TRADE/STREET NAME	WINDOW OF DETECTION
PSYCHIATRIC	CNS STIMULANTS / ADHD	Ritalinic Acid		50	Ritalin, Concerta	4d
		Atomoxetine		100	Strattera	7d
		Amphetamine		50	Adderall, Vyvanse	3d
		Methamphetamine		50	Desoxyn	3d
		Phentermine		100	Adipex-P, Suprenza	5d
	TRICYCLIC ANTIDEPRESSANTS	Clomipramine		50	Anafranil	14d
		Desipramine		100	Norpramin	2-11d
		Imipramine		50	Tofranil	5d
		Amitriptyline		50	Elavil	10d
		Nortriptyline		50	Pamelor	10d
		Doxepin	+	100	Sinequan	2-8d
		Protriptyline		100	Vivactil	27d
	Trimipramine		100	Surmontil	2d	
	BARBITURATES	Butalbital		400	Fioricet, Orbivan	7d
		Phenobarbital		400	Luminal, Solfoton	15d
		Secobarbital		400	Seconal	3d
	BENZODIAZEPINES/ ANXIOLYTICS	Alprazolam	+	50	Xanax	5d
		Clonazepam	+	50	Klonopin	5d
		Flunitrazepam	+	50	Rohypnol	5d
Flurazepam		+	50	Dalmane	2d	
Midazolam		+	50	Versed	2d	
Triazolam		+	50	Halcion	2d	
Diazepam		+	50	Valium	10d	
Chlordiazepoxide			50	Librium	5d	
Oxazepam			50	Serax	5d	
Estazolam			50	Prosom	4d	
Lorazepam			50	Ativan	7d	
Temazepam			50	Restoril	5d	
Buspirone			50	Buspar	4d	
HYPNOTICS		Zaleplon		50	Sonata	<2d
	Zolpidem		1	Ambien	4d	
PAIN / NEUROLOGY	OPIOIDS	Codeine		50	Tylenol #3	4d
		Morphine		50	MS Contin, Kadian	4d
		Hydrocodone		50	Lortab, Vicodin	4d
		Hydromorphone	+	50	Dilaudid, Exalgo	4d
		Oxycodone	+	50	Oxycontin, Percocet	4d
		Oxymorphone		50	Opana	4d
	SYNTHETIC OPIATES	Propoxyphene	+	50	Darvon, Darvocet	5d
		Buprenorphine	+	5	Suboxone, Butrans	7d
		Fentanyl	+	5	Actiq, Duragesic	3d
		Sufentanil		50	Sufenta	2d
		Tramadol	+	50	Ultram	4d
		Desomorphine		50	Krokodil	<1d
		Methadone	+	50	Methadone	5d
		Meperidine	+	50	Demerol, Pethidine	3d
Pentazocine		50	Talwin	2d		
Tapentadol	+	50	Nucynta	3d		

	CLASS	COMPOUND	METABOLINK™	CUTOFF (ng/mL)	TRADE/STREET NAME	WINDOW OF DETECTION	
PAIN / NEUROLOGY	OPIOID ANTAGONIST	Naloxone		5	Evzio, Narcan	<1d	
		Naltrexone		20	Revia, Vivitrol	4d	
	MUSCLE RELAXERS	Baclofen		100	Kemstro	3d	
		Cyclobenzaprine	+	50	Flexeril	8d	
		Carisoprodol		100	Soma	<1d	
		Meprobamate		100	Equanil	3d	
	ANTICONVULSANTS	Gabapentin		400	Neurontin	2d	
		Pregabalin		400	Lyrica	2d	
	NSAID's	4-Hydroxydiclofenac		100	Voltaren	<2d	
		Ibuprofen		400	Advil, Motrin	<1d	
		Naproxen		200	Aleve, Naprosyn	10d	
	ANESTHETIC	Ketamine	+	100	Ketamine	2d	
	OTC / HERB	OVER THE COUNTER	Acetaminophen		200	Tylenol	2d
			Pseudoephedrine		100	Sudafed	5d
Dextromethorphan			+	100	Robitussin, Delsym	2-5d	
PLANT ALKALOID		Mitragynine	+	50	Kratom	3d	
		Caffeine		100	Diurex, Vivarin	1-2d	
		Continie		200	(Nicotine metabolite)	4d	
ILLICITS	ILLICITS	Cocaine	+	50	Cocaine	3d	
		PCP		10	Angel Dust	8d	
		MDA		100	Sally	3d	
		MDEA		100	Eve	3d	
		MDMA		100	Ecstasy, Molly	3d	
		THC-COOH		25	Cannabis, Marinol	5-30d	
		6-Acetylmorphine		20	(Heroin metabolite)	4d	
	PSYCHEDELICS	DMT		100	Illicit psychadelic	2-5d	
		DOB		100	Illicit psychadelic	2-5d	
		LSD		10	Acid	<1d	
	BATH SALTS	Methylone		50	Illicit cathianone	4d	
		Ethylone		50	Illicit cathianone	4d	
		MDPV		100	Illicit cathianone	4d	
		Mephedrone		100	Illicit cathianone	4d	
		Naphyrone		100	Illicit cathianone	4d	
		Alpha-PVP		100	Flakka, Gravel	4d	
	SYNTHETIC CANNABINOIDS	JWH-018		50	K2/Spice	3d	
		JWH-073		50	K2/Spice	3d	
		JWH-081		50	K2/Spice	3d	
		JWH-122		50	K2/Spice	3d	
		JWH-250		50	K2/Spice	3d	
EIOH	ETHANOL METABOLITES	Ethyl Glucuronide		500	(Ethanol metabolite)	3d	
		Ethyl Sulfate		250	(Ethanol metabolite)	3d	
PSYCHIATRIC	SSRI/SNRI	Paroxetine		50	Paxil	4d	
		Sertaline		50	Zoloft	4d	
		Citalopram		50	Celexa, Lexapro	4d	
		Venlafaxine		50	Effexor	4d	
		Fluoxetine	+	50	Prozac	4d	
		Duloxetine		50	Cymbalta	4d	
		Norquetiapine		50	Seroquel	4d	
	ANTIPSYCHOTICS	Aripiprazole		50	Abilify	4d	
		Olanzapine	+	50	Zyprexa	4d	
		9-Hydroxyrisperidone		50	Risperidal	4d	
		Ziprasidone		50	Geodon	4d	
		Thioridazine		50	Thioridazine	4d	

+ Parent compound AND metabolite concentrations are analyzed and reported.